

# Agility® Now Runs on MS Windows® 10 OS!

## Upgrade to Security and Peace of Mind

Dynex® Technologies has made it possible to upgrade your Agility operating system to Windows 10 after Microsoft phased out security updates support for Windows 7 OS. While your Agility system will continue to function and produce results for your ELISA testing with Windows 7 OS, it will do so without continued software and security updates and at a greater risk for viruses and malware. The best way to stay secure is on Windows 10.

Dynex has created an upgrade kit package for upgrading to Windows 10 operating system that includes hardware, software, and transition your database with integrity for a smooth and seamless upgrade of your Agility system. Our network of highly trained service engineers will upgrade your system for continued security updates and peace of mind.

New Agility systems will now come with Windows 10 OS preinstalled.

## Protection and Satisfaction

- Periodic system updates from Microsoft
- Insurance against viruses, malware, and ransomware
- Safeguard your patient data

## Improved Performance

- New motherboard with 6<sup>th</sup> Generation Intel® Core™ i7 Processor delivers high performance, improved data integrity, and lower power consumption
- 16 GB RAM for a faster more responsive operating system
- 500 GB Solid State hard drive for large data storage capacity and faster data retrieval



## Agility Windows® 10 Upgrade Kit Ordering Information

<b>Part #</b>	<b>Product Description</b>	<b>List Price</b>
18000040	Agility Windows 10 upgrade kit and installation by service personnel	\$7,500